

HeartLine™

The Best of Tele-Cardiology

Aerotel's HeartLine electrocardiogram (ECG) remote monitoring system is winning the hearts of cardiologists and cardiac patients around the world. Patients conduct normal, mobile lives, while at the same time they and their attending physicians know that a professional response to any change in their heart condition is only a phone call away.

HeartLine is designed for various diagnostic and monitoring service applications. It comprises a complete range of personal, 1 to 12 lead ECG monitoring devices for various remote diagnostic and emergency service applications. The recorded ECG signals are transmitted via a variety of transmission media to a medical center where a PC-based ECG receiving station receives the transmitted ECG, for display, storage and analysis.

The HeartLine System Consists of:

- HeartView™ 12L - 12-Lead ECG Recorder/Transmitter, Acoustic/Digital
- HeartView™ P-12/8 Plus - 12/8-Lead ECG Personal Recorder/Transmitter, Acoustic
- HeartView™ P-12/8 BT - 12/8-Lead ECG Personal Recorder/Transmitter, Bluetooth®
- Heart 2005A™ - Single-Lead Transtelephonic ECG Loop Event Recorder/Transmitter
- Heart 2006™ - Dual-Lead Transtelephonic ECG Loop Event Recorder/Transmitter
- HeartSmart™ AF - Auto-Triggered Single Lead Transtelephonic ECG Loop Event Recorder/Transmitter
- HeartOne™ - Single Lead Transtelephonic ECG Post Event Recorder/Transmitter
- HeartLine Receiving Station (HRS™) - Cutting-Edge Technology for ECG Acquisition for Medical Monitoring Centers
- HRS-Net™ - Web-Based Remote Access to HRS Data
- HeartLine Remote Receiver (HRR™) - PC-Based ECG Acquisition and Transmission Application
- Voice Receiving Module (VRM™) - Automatic Answering & Receiving of ECG Transmissions



Heart 2006™/2005A™
Single/Dual lead transtelephonic ECG loop event recorder/transmitter

HeartLine


HeartView™ P-12/8 plus

Transtelephonic/Wireless
personal 12 lead ECG
recorder/transmitter



HeartView™ 12L

Transtelephonic/Wireless
professional 12 lead ECG
recorder/transmitter



Call Center




Internet



Web Access

HeartOne™

Single lead ECG,
post event
recorder/transmitter



HeartSmart™ AF

Single lead self triggered
transtelephonic ECG loop
event recorder/transmitter



Aerotel Medical Systems (1998) Ltd.

5 Hazoref St. Holon 58856, Israel

Tel +972 3 559 3222

Fax +972 3 559 6111

info@aerotel.com



www.aerotel.com